# STERRAD® INCUBATOR (58°C)

## REF: 21005 (100-120V~)

## 21006 (230V~)

**ASP** ADVANCED STERILIZATION PRODUCTS®

a *Johnson & Johnson* company

Division of Ethicon, Inc.

33 Technology Drive, Irvine, CA 92618-9824

Authorized EC Representative

Ethicon GmbH

Oststraße 1, D-22844 Norderstedt

<b>USA</b> Irvine, CA 92618	<b>I</b> 00040 Pomezia, Roma
<b>F</b> 92787 Issy-les-Moulineaux	<b>E</b> 28042 Campo de las Naciones, Madrid
<b>D</b> 22844 Norderstedt	<b>GB</b> Ascot, SL5 9EY
<b>CH</b> 8957 Spreitenbach	<b>CDN</b> Johnson & Johnson Medical Products Markham, ON, L3R 0T5
<b>NL</b> 3800 AD Amersfoort	<b>J</b> ジョーンソン・エンド・ジョーンソン 株式会社 〒101-0065 東京都千代田区西神田 3 丁目 5 番 2 号
<b>GR</b> 15125 Maroussi, Athens	<b>A</b> 1190 Wien
<b>S</b> 19184 Sollentuna	<b>BR</b> Rodovia Presidente Dutra, Km 154 S.J. Campos-S.P 12240-908
<b>B</b> 1700 Dilbeek	<b>P</b> 2745-555 Barcarena
<b>CZ</b> Na Radosti 399 155 25 Praha 5 - Zličín	<b>PL</b> Johnson & Johnson Poland Sp. z o.o. Szyszkowa 20, 20-285 Warszawa
<b>SK</b> Jakubovo nám. 13 811 09 Bratislava	<b>H</b> Johnson & Johnson Kft. H-2045 Törökbálint, Tó Park

English	2
Français	5
Deutsch	8
Italiano	11
Español	14
Nederlands	17
Português	20
Dansk	23
Suomi	26
Norsk	29
Svenska	32
Estonian	35
Ελληνικά	38
Japanese	41
Polski	44
Hungarian	47
Slovak	50
Czech	53

## INSTRUCTIONS FOR USE

### INTENDED USE:

The STERRAD® Incubator is designed to incubate up to 18 biological indicator vials at 58°C ±2°C. The open well design provides viewing of the vial contents while incubating. The STERRAD Incubator is designed for use with the Advanced Sterilization Products' (ASP) STERRAD CycleSure™ Biological Indicator. Use only as directed to ensure proper product performance.

### PRINCIPLES OF OPERATION:

The electronically controlled incubator is factory set at 58°C. Heating is achieved through a thermally conductive path from element to block to test tube to media. Heating control is achieved by means of a precision sensor and automatic electronic controller calibrated to the media temperature. Qualified technicians using a calibrated thermometer can perform field temperature adjustment, within the adjustable range.

### INSTRUCTIONS:

Place on dry horizontal surface with rubber feet down. Depress the rocker switch to the "on" position, which illuminates the switch to green. Wait at least 1 hour. The temperature should read 58°C ±2°C. The incubator is now ready for use. Refer to the STERRAD CycleSure Biological Indicator's instructions for incubation time.

### CALIBRATION:

Bring the incubator to 58°C. Place a calibrated thermometer in the thermometer port located on the left side of the top tier. Close the cover. Wait 15 minutes before reading the thermometer. No recalibration is needed if the thermometer reading is 58 ±2°C. The recalibration is performed with a screwdriver adjustment located at the center rear of the inner case. A small flat bladed screwdriver is needed to adjust the recalibration potentiometer. A clockwise quarter turn will decrease the temperature approximately 2°C. Wait at least 1 hour between adjustments.

### OPERATIONAL PRECAUTIONS AND LIMITATIONS:

1. Use only STERRAD CycleSure Biological Indicator to achieve specified incubator performance.
2. Locate the STERRAD Incubator in a static air area. Do not locate in the air stream of a heating or cooling vent.
3. Incubator should be plugged into a properly grounded electrical outlet.
4. Incubator must be used with the included AC cordset only.

5. Disconnect power before cleaning or servicing. Servicing should be done by a qualified technician only.
6. Risk of electric shock. Do not open.
7. Do not use in an explosive atmosphere.
8. Do not touch white block surface when incubator is turned on. Surface is hot.
9. Close cover during incubation period.
10. Replace fuses with those of the same type and rating. (Applicable to Model 21006 only.) Fuse size: 5 mm diameter x 20 mm length; Fuse rating: 1 AMP, 250 Volt.

**CLEANING:**

Unplug STERRAD Incubator. Wipe down with a lightly dampened cloth of sodium hypochlorite, ensuring that no liquid gets in the incubator. Dry with a clean cloth. Contact ASP before using any cleaning solution not recommended by the manufacturer.

**PERFORMANCE CHARACTERISTICS AND SPECIFICATIONS:**

MODEL	VOLTAGE	AMPERAGE	POWER	FREQUENCY	TEMP SETPOINT
21005	100-120 ±10% V~	0.4A	33-48W	50/60 Hz	58 ±2°C
21006	230 ±10% V~	0.2A	44W	50/60 Hz	58 ±2°C

**ENVIRONMENTAL CONDITIONS:**

## Operating:

15° - 30°C; 10 - 85% relative humidity at 25°C, non-condensing  
 Installation Category II (overvoltage) in accordance with IEC 664  
 Pollution Degree 2 in accordance with IEC 664  
 Tested to an altitude of 2,000 meters

## Storage:

-40° - 70°C; 10 - 85% relative humidity

**THERMOMETER:** 30° - 65°C, INDICATION ACCURACY ±0.5°C

Instructions For Reuniting Separated Fluid Column of Your EVER-SAFE™ Thermometer.

**IMPORTANT: Store Thermometers in an upright position only. Handle Instruments With Care. Wear Safety Glasses and Gloves Before Proceeding.**

**Heating Method:** Heat the thermometer bulb in an upright position away from your face in warm liquid, air, or over a soft flame sufficient to allow the liquid column to rise slowly until the separated portion of the column enters the

expansion chamber at the top of the thermometer. Note that over filling the expansion chamber will break the thermometer. Tap the thermometer gently on the surface of a large rubber stopper in an upright position allowing the gas separating the column to rise above the column. Allow the thermometer to cool slowly in an upright position.

**Cooling Method:** Prepare a solution of shaved ice and salt or CO<sub>2</sub> (Dry-Ice) and alcohol. Place the thermometer bulb only in the solution. Keep the thermometer upright. Allow the liquid column to retreat into the bulb, swing the thermometer (Bulb Down) in the arc releasing the entrapped gas permitting it to escape above the column. Allow the thermometer to warm slowly in an upright position.

**LISTING:** CUL (100-120V~), EN61010-1 (230V~), UL (100-120V~)

**SYMBOL DESCRIPTION:**

~ AC Alternating Current    I ON (Power Supply)    ⚠ Caution    🔥 Hot Surface  
 ⚡ Protective Earth    O OFF (Power Supply)    ⎓ Fuse

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